



ACCESSORIES
FEEL LIKE **ON CLOUD NINE**

BRUMABA

OPERATING TABLE SYSTEMS

BRUMABA is a second-generation, family-owned company that develops, manufactures, and globally distributes operating tables, operating chairs, and patient transporters and appropriate accessories for specialty departments in hospitals, clinics, and doctor's offices.

BRUMABA'S CORE COMPETENCE:
MOBILE, SPECIALIZED OPERATING TABLES.

BRUMABA is a professionally competent contact and adviser for all matters concerning ergonomics and positioning on operating tables. We offer our customers products for specialized departments such as eye and oral surgery, OMS, and plastic surgery, dermatology, as well as for general operations, both for in-patients and out-patients.

BRUMABA'S CORE BENEFIT:
OPERATING TABLES FOR OPERATING WITHOUT
TIRING AND OPTIMUM PATIENT POSITIONING.

This complete catalog of accessories contains appropriate accessories for expanding and customizing your entire specialty operating table system according to medical discipline.



POSITIONING | **HEAD**

Head plate, large, 400 x 300 mm	6
Head plate, small, 320 x 240 mm	6
Horseshoe-shaped headrest	6
Horseshoe-shaped headrest, small	7
Horseshoe-shaped headrest with long clamping piece, 110 mm	7
Silicone head restraining belt for fastening on the headrest	7
Hair transplant support	8
Neck cushion, viscoelastic, black	8
Head plate “rubens pillow” with ball-and-socket joint	8
Headrest, detachable, GENIUS (surcharge)	9
Headrest, detachable, FOCUS, PRIMUS, VARIUS (surcharge)	9
Radiolucent headrest, detachable (surcharge)	9
External mounting with head plate, radiolucent	10
Microsurgery wrist rest with integrated mounting bracket	10
Microsurgery wrist rest (GENIUS EYE)	11
Lower arm support for Ophthalmology	11
Tongue holder, extended, with flat plate rail	11

POSITIONING | **ARM/SHOULDER/TORSO/LEG**

Armrest OP with fastening belt	12
Arm support lateral with fastening belt	12
Mounting bracket, swiveling	12
Mounting bracket, swiveling (GENIUS)	12
Mounting bracket, swiveling (VENUS)	12
Hand table with mounting bracket SPECIAL “XL”	13
Hand table telescopic support, single	13
Hand surgery table with double telescopic support, with mounting bracket SPECIAL “XL”	13
Back section cushion with arm support, viscous	14
Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt	14
Back section extension for shoulder arthroscopy	14
Side guard, hinged	15
Lateral positioning support, deep	15
Lateral positioning support, dorsal (back)	15
Body belt, two-part, with 2 mounting gliders	16
Body belt , with 2 mounting gliders	16

POSITIONING | **LEG**

Leg holder, OP	17
Mounting adapter with 90° redirection for leg holder	17
Leg holder with 2 ball-and-socket joints	17
Joint support	18
Arthroscopic surgery leg support	18
Padded insert for arthroscopic surgery leg support	19
Standing plate (1 pair)	19

POSITIONING PADDING

Lordosis support padding	20
Half roll padding, small	20
Half roll padding, large	20
Stomach positioning padding, trapezoid	21
Mattress, viscoelastic	21
Mattress, viscoelastic, short, with gynecological cutout	21

OPERATING TABLE | **EXTENSIONS AND WIDENERS**

Back section extension, narrow (OP)	22
Back section extension, narrow (GENIUS)	22
Back section extension, wide (OP)	22
Back section extension, wide (GENIUS)	22
Back section extension, extra wide	23
Back section extension (GENIUS EYE / GENIUS EYE PRO)	23
Foot section extension, wide (OP)	23
Foot section extension, wide (GENIUS)	23
Side widener/foot section extension, narrow (GENIUS)	24
Side extension/foot section extension, narrow (OP)	24
Side extension/foot section extension, narrow (GENIUS)	24
Side extension, radiolucent	24
Lower leg plate, combined (VARIUS)	25
Lower leg plate, combined, extended (VARIUS)	25
Lower leg plates, extended (1 pair) (VARIUS)	25

INSTRUMENT TABLES

Monitor holder, swiveling, on column	26
Instrument table, mobile, 600 x 400 mm	26
Instrument table, mobile, 700 x 500 mm	26
Instrument table on the column, 600 x 400 mm	27

Instrument table on the column, 700 x 500 mm	27
Tray, Aluminium, for instrument table	27

ANESTHESIA, INFUSION STAND

Anesthetic frame, bendable, foldaway	28
Anesthetic frame, bendable, with oxygen supply	28
Anesthesiology tube holder	29
Anesthesiology tube holder with remote control holder	29
Anesthesiology tube holder for surgery, front end	29
Anesthesiology tube holder, front end (GENIUS)	29
Anesthetic frame, universal	30
Infusion stand with ball-and-socket joint	30
Infusion stand on the column, swiveling	31
Infusion stand - swivel arm, extended (360 mm)	31
Infusion stand on instrument table	31

CONTROL ELEMENTS

Hand-held remote control	32
Remote control holder on infusion stand	32
Remote control holder, flexible	33
Remote control holder, short	33
Magnetic plate for remote control	33
Foot control, type IV	34
Foot control, type IV without memory	34
Pedestal for foot control, type IV	35
Foot control, type IV fixed on undercarriage	35
Foot control, type IV EYE, without memory	35

MOUNTING | **BRACKETS AND , ADAPTERS**

Mounting bracket	36
Mounting bracket SPECIAL “XL”	36
Mounting glider with clamping	37
Mounting adapter, 25 x 10 mm	37
Mounting adapter, strenghtended, 290 mm	37

MOUNTING | **MOUNTING RAILS**

Mounting rail, 250 mm head section (GENIUS)	38
Mounting rail, 340 mm foot section (GENIUS)	38

Mounting rail, 250 mm, massiv seat section (GENIUS)	38
Mounting rail, 340 mm, solid back and foot section (GENIUS)	39
Mounting rail head section (GENIUS EYE)	39
Mounting rail set (VENUS)	39

TRANSPORT HANDLES

Transport handles, foldaway, head end (1 pair) (GENIUS)	40
Transport handles foldaway, head end, (1 pair) (FOCUS, PRIMUS, VARIUS)	40
Transport handles foldaway, foot end (1 pair)	40

HYGIENE/PRESERVATION/CARE

EasyCover, disposable surgery bed drape, unsterile	41
Paper towel holder	41
Paper towel holder, 550 mm (VENUS)	41
Protective film for foot section,visco (FOCUS, PRIMUS, GENIUS, VENUS)	42
Protective film for splited foot section, visco (VARIUS) (1 pair)	42
Service kit	43

MAGNIFYING LAMP

Magnifying lamp with swivel arm	43
Magnifying lamp with swivel arm (VENUS)	43

ACCESSORY ORDER

Wall mounting rail	44
Gliders for wall mounting rail / accessories trolley	44
Wall mount for hand surgery table	44
Accessory stand	44
Accessory stand, tray	45
Accessories trolley	45
Transport trolley for batteries	45

HEAD PLATE, LARGE, 400 x 300 mm



ORDER NUMBER: Z1.003


PADDING DIMENSIONS (L x W x H):
400 x 300 x 50 mm

CLAMPING PIECE HEIGHT: **110 mm**

- For head positioning, lateral positioning
- Abdominal positioning in conscious patients
- Ball joint for seamless inclination, height-adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HEAD PLATE, SMALL, 320 x 240 mm



ORDER NUMBER: Z1.109

PADDING DIMENSIONS (L x W x H):
320 x 240 x 50 mm

CLAMPING PIECE HEIGHT: **110 mm**

- For head positioning, lateral positioning
- Abdominal positioning in conscious patients
- Ball joint for seamless inclination, height-adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HORSESHOE-SHAPED HEADREST
(INCLUDED AS STANDARD)



ORDER NUMBER: Z1.043

PADDING DIMENSIONS (L x W x H):
250 x 240 x 110 mm

CLAMPING PIECE HEIGHT: **65 mm**

- For a defined and stable head position
- Ball-and-socket joint, height adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HORSESHOE-SHAPED HEADREST, SMALL



ORDER NUMBER: Z1.043A

PADDING DIMENSIONS (L x W x H):
200 x 210 x 80 mm

CLAMPING PIECE HEIGHT: **65 mm**

- For good access to the ears
- For head positioning in children
- Ball-and-socket hinge, height adjustable,

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HORSESHOE-SHAPED HEADREST WITH LONG CLAMPING PIECE, 110 mm



ORDER NUMBER: Z1.249

CLAMPING PIECE HEIGHT: **110 mm**

- For Bekhterev's disease patients
- Defined and stable head position
- Ball-and-socket hinge, height adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

SILICONE HEAD RESTRAINING BELT FOR FASTENING ON THE HEADREST



ORDER NUMBER: Z1.189

BELT SIZE (L x W):
900 x 50 mm

- Hygienic silicone belt for the fixation of the patient's head

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HAIR TRANSPLANT SUPPORT



ORDER NUMBER: Z1.032

PADDING DIMENSIONS (Ø x L):
90 x 230 mm

- For optimum access to the upper and rear head area
- Ball-and-socket hinge, height adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

NECK CUSHION, VISCOELASTIC, BLACK



ORDER NUMBER: Z1.250

PADDING DIMENSIONS (Ø x L):
200 x 200 mm
CLAMPING PIECE HEIGHT: **65 mm**

- For extra soft neck positioning
- Ball-and-socket hinge, height adjustable

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HEAD PLATE “RUBENS PILLOW” WITH BALL-AND-SOCKET JOINT



ORDER NUMBER: X.000184

PADDING DIMENSIONS (L x W x H):
47 x 25 x 10 mm
CLAMPING PIECE HEIGHT: **110 mm**

- For defined head positioning
- Armrest for the surgeon during eye surgery in 12 o'clock position
- Ball-and-socket hinge, height adjustable,

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

HEADREST, DETACHABLE GENIUS (SURCHARGE)



ORDER NUMBER: X.000257

- For a quick change of the head section
- Good access for the anesthetist, e.g. during knee or gynecological surgery

Suitable for:
» all BRUMABA Operating Tables except “VENUS”, “GENIUS EYE” and “GENIUS EYE PRO”

HEADREST, DETACHABLE FOCUS, PRIMUS, VARIUS (SURCHARGE)



ORDER NUMBER: X.000209

- For a quick change of the head section
- Good access for the anesthetist, e.g. during knee or gynecological surgery
- Solid stainless steel construction

Suitable for:
» all BRUMABA Operating Tables except “VENUS”, “GENIUS EYE” and “GENIUS EYE PRO”

RADIOLUCENT HEADREST, DETACHABLE (SURCHARGE)




ORDER NUMBER: X.000235

OPTIONAL: **PADDING FOR HEAD SECTION RADIOLUCENT (AS SHOWN)**
ORDER NUMBER: E.002553

- For a quick change of the head section
- Good for intra-operative X-ray, e.g. ENT, OMS

Suitable for:
» all BRUMABA Operating Tables except “VENUS” and “GENIUS EYE”



EXTERNAL MOUNTING WITH HEAD PLATE, RADIOLUCENT **ORDER NUMBER: Z1.194**



- DIMENSION (L x W):
330 x 250 mm
- For X-ray of the head or the whole spinal column in inverse position
 - Mounting at foot end
 - Detachable padding
 - Angle adjustment by cogging

Suitable for:
» BRUMABA Operating Table “VARIUS”

Max. head weight: 20 kg
Max. patient weight for inverse set-up 200 kg

MICROSURGERY WRIST REST WITH INTEGRATED MOUNTING BRACKET **ORDER NUMBER: Z1.240**



- DIMENSIONS:
INTERIOR Ø: **325 mm**
EXTERIOR Ø: **440 mm**
PADDING WIDTH: **40 mm**
- For precise operating and operating without tiring
 - Ideal for microscope surgery
 - Ball-and-socket hinge, height adjustable
 - Soft padding

Suitable for:
» all BRUMABA Operating Tables except „VENUS“, „GENIUS EYE“ and „GENIUS EYE PRO“

MICROSURGERY WRIST REST (GENIUS EYE)

ORDER NUMBER: Z1.257



- DIMENSIONS:
INTERIOR Ø: **230 mm**
EXTERIOR Ø: **315 mm**
- For precise operating and operating without tiring
 - Ideal for microscope surgery
 - Ball-and-socket hinge, height adjustable

Suitable for:
» BRUMABA Operating Tables “GENIUS EYE” and “GENIUS EYE PRO”

LOWER ARM SUPPORT FOR OPHTHALMOLOGY

ORDER NUMBER: X.000013



- Ergonomic lower armrest for the surgeon during eye surgery in 12 o'clock position
- For operating without tiring during ophthalmological procedures
- Soft padding

Suitable for:
» BRUMABA Operating Tables „GENIUS EYE“ and „GENIUS EYE PRO“

TONGUE HOLDER, EXTENDED, WITH FLAT PLATE RAIL

ORDER NUMBER: X.000218



- DIMENSIONS (W x H):
460 x 390 – 670 mm
- OPTIONAL: **TONGUE HOLDER, SHORT, WITH RAISED PLATE RIM (AS SHOWN)**
- ORDER NUMBER: 21.061**
- For ENT
 - Large individual range by way of two ball-and-socket hinges
 - Chain for tongue and retractor holder
 - instrument tray
 - Height-adjustable

Suitable for:
» all BRUMABA Operating Tables except “VENUS” and “GENIUS EYE”

ARMREST OP WITH FASTENING BELT



ORDER NUMBER: Z1.001

PADDING DIMENSIONS (L x W x H):
490 x 140 x 40 mm

RESTRAINING BELT:
650 x 50 mm

- For ergonomic arm positioning
- Eccentric lever operates two ball-and-socket hinges, for fast, seamless adjustment
- Also for prone position
- Integrated impact protection
- Detachable padding, Silicone restraining belt

Suitable for:
» all BRUMABA Operating Tables

ARM SUPPORT LATERAL WITH FASTENING BELT



ORDER NUMBER: Z1.002

PADDING DIMENSIONS (L x W x H):
240 x 140 x 40 mm

- For slender arm apposition
- Folds down by way of swivel joint
- Detachable padding, Silicone restraining belt

Suitable for:
» all BRUMABA Operating Tables

MOUNTING BRACKET, SWIVELING



ORDER NUMBER: Z1.033 (AS SHOWN)

MOUNTING BRACKET FOR "GENIUS"
ORDER NUMBER: Z1.281

MOUNTING BRACKET FOR "VENUS"
ORDER NUMBER: Z1.297

- For keeping the patient in the prone position

Suitable for:
» BRUMABA Operating Tables "FOCUS", "PRIMUS" and "VARIUS"

HAND SURGERY TABLE WITH MOUNTING BRACKET SPECIAL "XL"



ORDER NUMBER: Z1.037

PADDING DIMENSION (L x W):
660 x 400 mm

- For ergonomic positioning of the patient's arm
- Ball-and-socket hinge
- Radiolucent
- Slim design for maximum legroom

Suitable for:
» all BRUMABA Operating Tables except "GENIUS EYE", "GENIUS EYE PRO" and "VENUS"

HAND TABLE - TELESCOPIC SUPPORT, SINGLE



ORDER NUMBER: Z1.089

HEIGHT ADJUSTABLE FROM:
570 - 1030 mm

- For greater stability requirements at the hand table **Order number: Z1.037**
- Height-adjustable

Suitable for:
» all BRUMABA Operating Tables except "GENIUS EYE", "GENIUS EYE PRO" and "VENUS"

HAND SURGERY TABLE WITH DOUBLE TELESCOPIC SUPPORT, WITH MOUNTING BRACKET SPECIAL "XL"



ORDER NUMBER: Z1.036

PADDING DIMENSION (L x W):
800 x 500 mm

HEIGHT ADJUSTABLE FROM:
690 - 1110 mm

- For hand surgery with greater space requirements
- Radiolucent
- Ball-and-socket hinge
- Two height-adjustable telescopic supports for extra steadiness

Suitable for:
» all BRUMABA Operating Tables except "GENIUS EYE", "GENIUS EYE PRO" and "VENUS"

BACK SECTION CUSHION WITH ARM SUPPORT, VISCOUS

ORDER NUMBER: Z1.259



PADDING WIDTH:
720 mm


- For simple arm positioning

Suitable for:

» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“

HORSESHOE-SHAPED HEADREST FOR SHOULDER ARTHROSCOPY WITH HEAD RESTRAINING BELT

ORDER NUMBER: Z1.154



PADDING DIMENSIONS (L x W x H):
200 x 210 x 80 mm

CLAMPING PIECE HEIGHT: **110 mm**

- For laterally offset positioning of the head
- For free access to the left and right shoulder area
- In combination with the head cradle, the horsehoe-shaped headrest adapts automatically during adjustment
- Ball-and-socket hinge, height adjustable, silicone restraining belt

Suitable for:

» BRUMABA Operating Tables „FOCUS“, „PRIMUS“ and „VARIUS“

BACK SECTION EXTENSION FOR SHOULDER ARTHROSCOPY

ORDER NUMBER: Z1.021



PADDING DIMENSIONS (L x W x H):
270 x 150 x 60 mm

- For shoulder operations left and right
- Easy assembly

Suitable for:

» BRUMABA Operating Tables „FOCUS“, „PRIMUS“ and „VARIUS“

SIDE GUARD, HINGED

ORDER NUMBER: Z1.187



DIMENSIONS (W x H):
350 x 310 mm

- For additional safety during transport of the patient
- Folds down by way of swivel joint

Suitable for:

» all BRUMABA Operating Tables except „VENUS“

LATERAL POSITIONING SUPPORT, DEEP

ORDER NUMBER: Z1.110



PADDING DIMENSIONS (L x W x H):
180 x 300 x 50 mm

CLAMPING PIECE HEIGHT: **120 mm**

OPTIONAL:
LATERAL POSITIONING SUPPORT
ORDER NUMBER: Z1.054
CLAMPING PIECE HEIGHT: **65 mm**

- For stable lateral positioning of the patient
- Ball-and-socket hinge for seamless adjustment

Suitable for:

» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“

LATERAL POSITIONING SUPPORT, DORSAL (BACK)

ORDER NUMBER: Z1.152



PADDING DIMENSIONS (L x W x H):
240 x 140 x 40 mm

HEIGHT-ADJUSTABILITY **120 mm**

- For slender support of the patient in lateral position
- Two joints for individual adaptation
- Height-adjustable

Suitable for:

» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“



BODY BELT, TWO-PART, WITH 2 MOUNTING GLIDERS

ORDER NUMBER: Z1.241



DIMENSION (L x W):
1475 x 80 mm

- For securing the patient
- Hygienic silicon material
- Reduced pressure load

Suitable for:
» all BRUMABA Operating Tables

BODY BELT, WITH 2 MOUNTING GLIDERS

**ORDER NUMBER: Z1.112
E.001061 | GENIUS EYE**



DIMENSION (L x W):
1440 x 80 mm

- For securing the patient
- Hygienic silicon material
- Reduced pressure load

Suitable for:
» all BRUMABA Operating Tables

LEG HOLDER OP



ORDER NUMBER: Z1.087

PADDING DIMENSIONS (L x W x H):
220 x 330 X 90 mm
RESTRAINING BELT DIMENSIONS (L x B):
650 x 50 mm

HEIGHT-ADJUSTABILITY: **130 mm**

- For ergonomic leg positioning
- Three-dimensional mobility
- Easy to handle by just one person
- Rocker function for automatic adaptation of the leg band against positional damage
- Ball-and-socket joints, height adjustable, silicon restraining belt

Suitable for:
» all BRUMABA OP Tables except „GENIUS EYE“ and „GENIUS EYE PRO“

MOUNTING ADAPTER WITH 90° DEFLECTION FOR SURGICAL LEG HOLDER



- For deflection of the surgical leg holder
Order number: Z1.087
for side positioning

Suitable for:
» all BRUMABA Operating Tables except “GENIUS EYE” and “GENIUS EYE PRO”

LEG HOLDER WITH 2 BALL-AND-SOCKET JOINTS




ORDER NUMBER: X.000126

PADDING DIMENSIONS (L x W x H):
220 x 330 X 90 mm

HEIGHT-ADJUSTABILITY **170 mm**

- For ergonomic leg positioning
- Three-dimensional mobility
- Easy to handle by just one person
- Suitable for side positioning
- Two ball-and-socket joints
- Height-adjustable
- Silicone restraining belt

Suitable for:
» all BRUMABA OP Tables except „GENIUS EYE“ and „GENIUS EYE PRO“



JOINT SUPPORT

ORDER NUMBER: Z1.028



PADDING DIMENSIONS (Ø x L):
90 x 230 mm
HEIGHT ADJUSTABLE FROM:
380 – 480 mm

- For surgery on the foot or shin
- Ball-and-socket joint
- Height-adjustable

Suitable for:
» all BRUMABA Operating Tables

ARTHROSCOPIC SURGERY LEG SUPPORT


ORDER NUMBER: Z1.016



MAX. Ø THIGH:
250 mm
INCL. MOUNTING BRACKET XL
ORDER NUMBER: Z1.088
INCL. MOUNTING BRACKET
ORDER NUMBER: Z1.047


- for knee arthroscopic surgery
- Suitable for left and right knee
- Simple thigh adaptation
- Also for cruciate ligament surgery
- Can be rotated 360°
- Height-adjustable

Suitable for:
» BRUMABA Operating Table „VARIUS“



**PADDED INSERT FOR ARTHROSCOPIC SURGERY
LEG SUPPORT**

ORDER NUMBER: Z1.155




PADDING THICKNESS: **25 mm**

- For knee arthroscopic surgery
without tourniquet

Suitable for:
» BRUMABA Operating Table „VARIUS“

STANDING PLATE (1 PAIR)

ORDER NUMBER: Z1.057



PADDING DIMENSIONS (B x L):
180 x 300 mm

- For procedures in the reverse
Trendelenburg position, e.g. phlebology

Suitable for:
» BRUMABA Operating Tables “VARIUS”
and “VENUS”

FOR “PRIMUS”
ORDER NUMBER: x.000359

POSITIONING PADDINGS

LORDOSIS SUPPORT PADDING



ORDER NUMBER: Z1.140

PADDING DIMENSIONS (L x W x H):
400 x 210 x 60 mm

- Soft material and hygienic, seamless cover
- For the support of the lordosis

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, SMALL



ORDER NUMBER: Z1.100

PADDING DIMENSIONS (L x W x H):
400 x 140 x 70 mm

- Soft material and hygienic, seamless cover
- Comfortable shape for various purposes

Suitable for:
» all BRUMABA Operating Tables

HALF ROLL PADDING, LARGE



ORDER NUMBER: Z1.101

PADDING DIMENSIONS (L x W x H):
555 x 210 x 120 mm

- Soft material and hygienic, seamless cover
- Various purposes

Suitable for:
» all BRUMABA Operating Tables

STOMACH POSITIONING PADDING, TRAPEZOIDAL (1 PIECE)



ORDER NUMBER: Z1.098

PADDING DIMENSION (L x W):
220 x 300 mm

MAX. DEPTH: **110 mm**

- Soft material and hygienic, seamless cover
- For ergonomic stomach positioning to relieve the abdominal area

Suitable for:
» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC



ORDER NUMBER: Z1.171

PADDING DIMENSIONS (L x W x H):
1660 x 570 x 40 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:
» all BRUMABA Operating Tables

MATTRESS, VISCOELASTIC, SHORT, WITH GYNECOLOGICAL CUTOUT




ORDER NUMBER: Z1.213

PADDING DIMENSIONS (L x W x H):
950 x 550 mm

- For long operations to prevent bedsores and for extra comfort
- Anti-static

Suitable for:
» BRUMABA Operating Table „VARIUS“



BACK SECTION EXTENSION, NARROW

ORDER NUMBER: Z1.006 | OP
Z1.183 | GENIUS



PADDING DIMENSION (H x W):
100 x 420 mm

- For correct positioning of the shoulders in large patients
- Included as standard with FOCUS, PRIMUS, VARIUS and GENIUS

Suitable for:
» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“

BACK SECTION EXTENSION, WIDE

ORDER NUMBER: Z1.005 | OP
Z1.184 | GENIUS



PADDING DIMENSION (H x W):
150 x 420 mm

- For correct positioning of the shoulders in large patients

Suitable for:
» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“

BACK SECTION EXTENSION, EXTRA WIDE

ORDER NUMBER: Z1.212

PADDING DIMENSION (H x W):
150 x 550 mm

- For correct positioning of the shoulders in large patients
- Extra wide

Suitable for:
» all BRUMABA Operating Tables except „GENIUS EYE“, „GENIUS EYE PRO“ and „VENUS“



BACK SECTION EXTENSION
(GENIUS EYE / GENIUS EYE PRO)

ORDER NUMBER: Z1.279

PADDING DIMENSION (H x W):
150 x 530 mm

- For correct positioning of the shoulders in large patients

Suitable for:
» BRUMABA Operating Tables „GENIUS EYE“ and „GENIUS EYE PRO“



FOOT SECTION EXTENSION, WIDE

ORDER NUMBER: Z1.268 | OP
Z1.286 | GENIUS*




PADDING DIMENSION (H x W):
150 x 420 mm

- For particularly large patients
- With edge protection plate
- Additional mounting bracket required

Suitable for:
» all BRUMABA Operating Tables except „VARIUS“

* INCL. MOUNTING ADAPTER



SIDE EXTENSION / FOOT SECTION EXTENSION, NARROW

ORDER NUMBER: Z1.291 | OP Z1.263 | GENIUS*



DIMENSIONS (H x W):
100 x 420 mm

OPTIONAL:
WITH MOUNTING RAIL
ORDER NUMBER: X.000644

■ For particularly large patients or widening of the lying surface

Suitable for:
» all BRUMABA Operating Tables

* INCL. MOUNTING ADAPTER

SIDE EXTENSION, RADIOLUCENT

ORDER NUMBER: Z1.055



DIMENSION (L x W x H):
410 x 225 X 55 mm

■ For evacuation of the extremities during operations and X-rays

■ Additional mounting bracket required

Suitable for:
» all BRUMABA OP Tables except „GENIUS EYE“ and „GENIUS EYE PRO“

LOWER LEG PLATE, COMBINED (VARIUS)

ORDER NUMBER: Z1.170



PADDING DIMENSIONS | VSKO(L x W x H):
500 x 555 x 60 mm

■ For use during inverse positioning e.g. with external mounting
Order number: Z1.194 for a homogeneous X-ray image

■ Connecting plate for the two-part foot section

Suitable for:
» BRUMABA Operating Table „VARIUS“

LOWER LEG PLATE, COMBINED, EXTENDED (VARIUS)

ORDER NUMBER: Z1.169



PADDING DIMENSIONS | VSKO (L x W x H):
800 x 555 x 60 mm

■ For use during inverse positioning

■ Extension of the X-ray range

■ For a homogeneous X-ray image

Suitable for:
» BRUMABA Operating Table „VARIUS“

LOWER LEG PLATES, EXTENDED, 1 PAIR (VARIUS)

ORDER NUMBER: Z1.020




PADDING DIMENSIONS (L x W x H):
800 x 270 x 50 mm

■ Extension of the leg plates e.g. after lithotomy position or knee surgery

Suitable for:
» BRUMABA Operating Table „VARIUS“

MONITOR HOLDER/ INSTRUMENT TABLE



MONITOR HOLDER, SWIFELING, ON COLUMN

ORDER NUMBER: X.000097



PLATE SIZE: 400 x 360 mm

- Not height-adjustable

Suitable for:
» all BRUMABA Operating Tables

INSTRUMENT TABLE, MOBILE

ORDER NUMBER: Z1.172




PLATE SIZE: 600 x 400 mm

OPTIONAL:
ORDER NUMBER: Z1.282
PLATE SIZE: 700 x 500 mm
HEIGHT ADJUSTABLE FROM:
560 - 1100 mm

- For the storage of instruments within easy reach
- Mobile with hygienic, smooth-running double castors
- One-handed height adjustment via gas spring
- Brake against rolling away

INSTRUMENT TABLES



INSTRUMENT TABLE ON COLUMN

ORDER NUMBER: Z1.275



PLATE SIZE: 600 x 400 mm

OPTIONAL:

ORDER NUMBER: Z1.278
PLATE SIZE: 700 x 500 mm

HEIGHT-ADJUSTABLE ABOVE
TABLE LEVEL: 300 - 540 mm

- For the storage of instruments within ergonomic reach above the patient
- One-handed adjustment via gas spring
- Plate 360° rotatable, arm 180° tiltable
- Left and right mounting possible

Suitable for:
» all BRUMABA Operating Tables

TRAY, ALUMINIUM, FOR INSTRUMENT TABLE

ORDER NUMBER: Z1.245




PLATE SIZE: 565 x 365 mm

- For efficient surgery organization
- Removable plate for the preparation of your instruments
- For instrument table
Order number: Z1.172 / Z1.275

Suitable for:
» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STAND



ANESTHETIC FRAME, BENDABLE, FOLDABLE

ORDER NUMBER: Z1.273



LENGTH: 960 mm

- For the supply of oxygen to the patient's head
- External hygienic tube guide on clips for 6 mm and 18 mm tubes (recommendation)
- Spacer for surgical drape
- Flexible gooseneck for individual adjustment
- Folds down by way of swivel joint

Suitable for:
» all BRUMABA Operating Tables

ANESTHETIC FRAME, BENDABLE, WITH OXYGEN SUPPLY

ORDER NUMBER: Z1.248



LENGTH: 1010 mm


- For the supply of oxygen to the patient's head
- Oxygen tube connection at gooseneck with outlet at head
- Spacer for surgical drape
- Flexible gooseneck for individual adjustment

Suitable for:
» all BRUMABA Operating Tables

ANESTHESIA / INFUSION STAND

ANESTHESIOLOGY TUBE HOLDER

ORDER NUMBER: Z1.049



LENGTH: 300 mm + CROWN

- For holding all standard tube diameters (similar to the Ulmer wheel)
- Flexible gooseneck for individual adjustment,

Suitable for:
» all BRUMABA Operating Tables

ANESTHESIOLOGY TUBE HOLDER WITH REMOTE CONTROL HOLDER

ORDER NUMBER: X.000322




- For anesthesiology tube holder
Order number: Z1.049

Suitable for:
» all BRUMABA Operating Tables


ANESTHESIOLOGY TUBE HOLDER, FRONT END

ORDER NUMBER: X.000256 | OP
X.000172 | GENIUS



- For front end tube guide
- Nasal anesthesia
- For holding all standard tube diameters (similar to the Ulmer wheel)

Suitable for:
» BRUMABA Operating Tables "FOCUS",
"PRIMUS" and "VARIUS"



ANESTHETIC FRAME, UNIVERSAL

ORDER NUMBER: Z1.015



WIDTH: 680 – 1130 mm
HEIGHT-ADJUSTABLE: 270 mm

- For all standard covering systems in surgery
- Length and height adjustable
- Ball-and-socket hinge
- Hook for infusion bottle

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND WITH BALL-AND-SOCKET JOINT

ORDER NUMBER: Z1.038



HEIGHT: 1110 – 1680 mm
WITH 4 HOOKS
MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on rail system
- One-handed height adjustment
- Ball-and-socket hinge
- 4 bottle holders

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND ON COLUMN, SWIVELING

ORDER NUMBER: Z1.181



HEIGHT: 1110 – 1680 mm
WITH 4 HOOKS
MAX. PERMISSIBLE LOAD:
2 KG PER HOOK

- Mounting on the column
- Independent of the angle of inclination setting of the patient surface
- One-handed height adjustment
- Tilttable
- 4 bottle holders

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND - SWIVEL ARM, EXTENDED (360 mm)

ORDER NUMBER: X.000214



OPTIONAL TO ORDER NUMBER: Z1.181
LENGTH: 360 mm

- For a larger tilting range of the infusion stand
ORDER NUMBER: Z1.181

Suitable for:
» all BRUMABA Operating Tables

INFUSION STAND ON INSTRUMENT TABLE


ORDER NUMBER: Z1.280



- Mounting on instrument table on column
ORDER NUMBER: Z1.275 / Z1.278
- One-handed height adjustment
- 2 bottle holders


Suitable for:
» all BRUMABA Operating Tables

CONTROLS



HAND-HELD REMOTE CONTROL

ORDER NUMBER: Z1.118



- Included as standard with delivery of the operating table
- Magnetic back
- Carry out all individual settings
- Retrieval of all 8 memory positions

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

REMOTE CONTROL HOLDER ON INFUSION STAND

ORDER NUMBER: Z1.243



- For attaching the hand-held remote control to the infusion stand

Suitable for:
» all BRUMABA Operating Tables

CONTROLS

REMOTE CONTROL HOLDER, FLEXIBLE

ORDER NUMBER: Z1.149



GOOSENECK LENGTH : 300 mm

- Flexible gooseneck

Suitable for:
» all BRUMABA Operating Tables

REMOTE CONTROL HOLDER, SHORT

ORDER NUMBER: Z1.026



- Mounting on rail system underneath patient surface

Suitable for:
» all BRUMABA Operating Tables

MAGNETIC PLATE FOR REMOTE CONTROL


ORDER NUMBER: Z1.045



- Included as standard with delivery of every operating table
- Attachment by way of rear adhesive surface

Suitable for:
» all BRUMABA Operating Tables

CONTROLS



FOOT CONTROL, TYPE IV



ORDER NUMBER: Z1.185.01
ORDER NUMBER: Z1.303 | VENUS

WEIGHT: 4.25 kg

- Carry out all individual settings
- Retrieval of all 8 memory positions
- Lifting bar for lifting or moving the foot control with the foot
- Can be operated in sterile condition

Suitable for:
» all BRUMABA Operating Tables except „VENUS“

FOOT CONTROL, TYPE IV WITHOUT MEMORY



ORDER NUMBER: Z1.271.01

WEIGHT: 4.25 kg

- Carry out all individual settings
- Lifting bar for lifting or moving the foot control with the foot
- No retrievable memory positions

Suitable for:
» all BRUMABA Operating Tables

PEDESTAL FOR FOOT CONTROL, TYPE IV




ORDER NUMBER: Z1.274

- For positioning the foot control on the undercarriage 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:
» all BRUMABA Operating Tables

FOOT CONTROL, TYPE IV FIXED ON UNDERCARRIAGE



ORDER NUMBER: Z1.225.01

WEIGHT: 4.25 kg

- Carry out all individual settings
- Retrieval of all 8 memory positions
- Foot control fixed on undercarriage 4-EASY-DRIVE, 5-WHEEL-DRIVE

Suitable for:
» all BRUMABA Operating Tables

FOOT CONTROL, TYPE IV EYE, WITHOUT MEMORY




ORDER NUMBER: Z1.226.01

WEIGHT: 4.25 kg

- Adjustment of head section, back section and height possible
- Foot control fixed on undercarriage
- No retrievable memory positions

Suitable for:
» all BRUMABA Operating Tables



MOUNTING BRACKET

ORDER NUMBER: **Z1.004**



For fastening of BRUMABA accessories

Suitable for:
» all BRUMABA Operating Tables

MOUNTING BRACKET SPECIAL “XL”

ORDER NUMBER: **Z1.088**



For fastening hand tables and arthroscopy leg holders

Suitable for:
» all BRUMABA Operating Tables

MOUNTING GLIDER WITH CLAMPING

ORDER NUMBER: **Z1.093**



For the permanent fastening of accessories

Suitable for:
» all BRUMABA Operating Tables

MOUNTING ADAPTER, 25 x 10 mm

ORDER NUMBER: **Z1.047**



DIMENSIONS (L x H x B):
200 x 25 x 10 mm
OPTIONAL: MOUNTING ADAPTER,
30 x 10 mm (SPECIAL)
ORDER NUMBER: Z1.047A
DIMENSIONS (L x H x B):
200 x 30 x 10 mm

For fastening accessories of other manufacturers (standard rails)

Suitable for:
» all BRUMABA Operating Tables

MOUNTING ADAPTER, STRENGTHENED, 290 mm

ORDER NUMBER: **X.000244**




DIMENSIONS (L x H x B):
290 x 25 x 10 mm

For fastening accessories of other manufacturers (standard rails)

Suitable for:
» all BRUMABA Operating Tables

MOUNTING RAIL, 250 mm HEAD SECTION (GENIUS)

ORDER NUMBER: Z1.151



■ For head section left

Suitable for:
» BRUMABA Operating Table “GENIUS“

MOUNTING RAIL, 340 mm FOOT SECTION (GENIUS)

ORDER NUMBER: Z1.079




■ For foot section left or right

Suitable for:
» BRUMABA Operating Table „GENIUS“

MOUNTING RAIL, 250 mm MASSIVE SEAT SECTION (GENIUS)

ORDER NUMBER: X.000087



■ For seat section, to attach leg holders
ORDER number: Z1.087

Suitable for:
» BRUMABA Operating Table “GENIUS“

MOUNTING RAIL, 340 mm SOLID
BACK AND FOOT SECTION (GENIUS)

ORDER NUMBER: X.000036




■ For back and foot section

Suitable for:
» BRUMABA Operating Table „GENIUS“

MOUNTING RAIL, 125 mm HEAD SECTION (GENIUS EYE)

ORDER NUMBER: Z1.262




■ E.g. for mounting the
microsurgery ring GENIUS EYE
ORDER NUMBER: Z1.257

Suitable for:
» BRUMABA Operating Tables “GENIUS EYE”
and “GENIUS EYE PRO”

MOUNTING RAILS SET (VENUS)

ORDER NUMBER: Z1.304




■ For back, seat and foot section

Suitable for:
» BRUMABA Operating Table “VENUS”

TRANSPORT HANDLES

TRANSPORT HANDLES, FOLDAWAY, HEAD END, GENIUS (1 PAIR)



ORDER NUMBER: Z1.210

- For the back-friendly moving of the operating table in chair position
- Foldaway for legroom

Suitable for:
» BRUMABA Operating Table “GENIUS”,
» with “GENIUS EYE” and “GENIUS EYE PRO”
included as standard

TRANSPORT HANDLES, FOLDAWAY, HEAD END, FOCUS, PRIMUS, VARIUS (1 PAIR)



ORDER NUMBER: Z1.211

- For the back-friendly moving of the operating table in chair position

Suitable for:
» BRUMABA Operating Tables “FOCUS”,
“PRIMUS” and “VARIUS”

TRANSPORT HANDLES, FOLDAWAY, FOOT END (1 PAIR)




ORDER NUMBER: Z1.164

- For the back-friendly moving of the operating table from the foot end

Suitable for:
» all BRUMABA Operating Tables

HYGIENE/PRESERVATION/CARE



EASYCOVER, DISPOSABLE SURGERY BED DRAPE, UNSTERILE

ORDER NUMBER: Z1.138



- DIMENSION (L x W):
2000 x 720 mm
- PACKAGING UNIT: **150 PIECES**
COLOR: **WHITE**
- Watertight disposable surgery bed drape for maximum hygiene
 - Head pocket prevents cover from slipping
 - Paper coating for comfort while lying

Suitable for:
» all BRUMABA Operating Tables

PAPER TOWEL HOLDER

ORDER NUMBER: Z1.051




- DIMENSIONS (W):
560mm
- USE OF PAPER TOWELS UP TO
MAX. WIDTH OF 560 mm
- For attaching paper towels under the foot section

Suitable for:
» all BRUMABA Operating Tables except
„VENUS“

VENUS ORDER NUMBER: Z1.294

DIMENSIONS (W):
550mm

PROTECTIVE FILM FOR FOOT SECTION, VISCO



ORDER NUMBER: Z1.223 | FOCUS/PRIMUS

- For the mechanical protection of the foot section padding e.g. outdoor footwear
- Easy disassembly for cleaning

Suitable for:
» BRUMABA Operating Tables „FOCUS“, „PRIMUS“

Optional:
ORDER NUMBER: Z1.224 | GENIUS
Z1.299 | VENUS

PROTECTIVE FILM FOR SPLITED FOOT SECTION, VISCO (VARIUS) (1 PAIR)



ORDER NUMBER: Z1.148A

- For the mechanical protection of the foot section padding e.g. outdoor footwear
- Easy disassembly for cleaning

Suitable for:
» BRUMABA Operating Table „VARIUS“

SERVICE KIT



ORDER NUMBER: Z1.305

- Everything to care for and clean your BRUMABA operating table
- Matting sponge, sponge, microfiber cloth, lubricant, cold welding paste, upholstery care, stainless steel care

Suitable for:
» all BRUMABA products

MAGNIFYING LAMP WITH SWIVEL ARM



ORDER NUMBER: Z1.217

LENS DIMENSION Ø: 160 mm
LENS DIOPTR 3.5
COLOR RENDITION: 6500 K
LED ENERGY EFFICIENCY: A+

Suitable for:
» all BRUMABA Operating Tables except “GENIUS EYE” and “VENUS”

Optional:
ORDER NUMBER: Z1.300 | VENUS

WALL MOUNTING RAIL



ORDER NUMBER: Z1.162

LENGTH: 100 mm

- For fastening accessories to the wall
- 5 gliders included

Suitable for:
» all BRUMABA Operating Table accessories

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:
ORDER NUMBER: Z1.161

WALL MOUNT FOR HAND SURGERY TABLE



ORDER NUMBER: Z1.190

- For space-saving storage of the hand table

Suitable for:
» all BRUMABA Operating Tables

ACCESSORY STAND



ORDER NUMBER: Z1.117

HEIGHT: 1080 mm
Ø BASE: 620 mm

- For storage of up to 8 accessories
- Mobile

Suitable for:
» all BRUMABA Operating Tables

Optional accessory:
Tray ORDER NUMBER: Z1.153

ACCESSORY STAND TRAY



ORDER NUMBER: Z1.153

Ø TRAY: 620 mm

- Additional storage option
for the accessory stand
ORDER NUMBER: Z1.117

Suitable for:
» all BRUMABA Operating Tables

ACCESSORIES TROLLEY



ORDER NUMBER: Z1.060

DIMENSIONS (L x W x H):
730 x 550 x 1000 mm
RAIL LENGTH: 725 mm

- For the storage of accessories
- Mobile
- 12 gliders included
- Extension for hand table and arthroscopic surgery leg support

Suitable for:
» all BRUMABA Operating Tables

GLIDERS FOR WALL MOUNTING /
ACCESSORIES TROLLEY:
ORDER NUMBER: Z1.161

TRANSPORT TROLLEY FOR BATTERIES



ORDER NUMBER: Z1.175

DIMENSIONS (L x W x H):
300 x 450 x 800 mm

- For easy transportation of two battery packs
- Holder for the charger
ORDER NUMBER: Z1.134

Suitable for:
» BRUMABA Operating Tables „FOCUS“, „PRIMUS“ and „VARIUS“

A

Accessory stand _____ 44

Accessory stand, tray _____ 45

Accessories trolley _____ 45

Anesthetic frame, bendable, foldaway _____ 28

Anesthetic frame, bendable, with oxygen supply _____ 28

Anesthetic frame, universal _____ 30

Anesthesiology tube holder _____ 29

Anesthesiology tube holder with remote control holder _____ 29

Anesthesiology tube holder for surgery, front end _____ 29

Anesthesiology tube holder, front end (GENIUS) _____ 29

Armrest OP with fastening belt _____ 12

Arm support lateral with fastening belt _____ 12

Arthroscopic surgery leg support _____ 18

B

Back section cushion with arm support, viscous _____ 14

Back section extension (GENIUS EYE / GENIUS EYE PRO) _____ 23

Back section extension for shoulder arthroscopy _____ 14

Back section extension, narrow (OP) _____ 22

Back section extension, narrow (GENIUS) _____ 22

Back section extension, wide (OP) _____ 22

Back section extension, wide (GENIUS) _____ 22

Back section extension, extra wide _____ 23

Body belt, two-part, with 2 mounting gliders _____ 16

Body belt , with 2 mounting gliders _____ 16

E

EasyCover, disposable surgery bed drape, unsterile _____ 41

F

Foot control, type IV _____ 34

Foot control, type IV fixed on undercarriage _____ 35

Foot control, type IV EYE, without memory _____ 35

Foot control, type IV without memory _____ 34

Foot section extension, wide (GENIUS) _____ 23

Foot section extension, wide (OP) _____ 23

G

Gliders for wall mounting rail / accessories trolley _____ 44

H

Hair transplant support _____ 8

Half roll padding, small _____ 20

Half roll padding, large _____ 20

Hand-held remote control _____ 32

Hand table with mounting bracket SPECIAL “XL” _____ 13

Hand table telescopic support, single _____ 13

Hand surgery table with double telescopic support, with mounting bracket SPECIAL “XL” _____ 13

Head plate, large, 400 x 300 mm _____ 6

Head plate, small, 320 x 240 mm _____ 6

Head plate “rubens pillow” with ball-and-socket joint _____ 8

Headrest, detachable, GENIUS (surcharge) _____ 9

Headrest, detachable, FOCUS, PRIMUS, VARIUS (surcharge) _____ 9

Horseshoe-shaped headrest _____ 6

Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt _____ 14

Horseshoe-shaped headrest for shoulder arthroscopy with head restraining belt _____ 14

Horseshoe-shaped headrest, small _____ 7

Horseshoe-shaped headrest with long clamping piece, 110 mm _____ 7

I

Infusion stand - swivel arm, extended (360 mm) _____ 31

Infusion stand on the column, swiveling _____ 31

Infusion stand with ball-and-socket joint _____ 30

Infusion stand on instrument table _____ 31

Instrument table, mobile, 600 x 400 mm _____ 26

Instrument table, mobile, 700 x 500 mm _____ 26

Instrument table on the column, 600 x 400 mm _____ 27

Instrument table on the column, 700 x 500 mm _____ 27

J

Joint support _____ 18

L

Lateral positioning support, deep _____ 15

Lateral positioning support, dorsal (back) _____ 15

Leg holder with 2 ball-and-socket joints _____ 17

Leg holder, OP _____ 17

Lower arm support for Ophthalmology _____ 11

Lower leg plate, combined (VARIUS) _____ 25

Lower leg plate, combined, extended (VARIUS) _____ 25

Lower leg plates, extended (1 pair) (VARIUS) _____ 25

Lordosis support padding _____ 20

M

Magnetic plate for remote control _____ 33

Magnifying lamp with swivel arm _____ 43

Magnifying lamp with swivel arm (VENUS) _____ 43

Mattress, viscoelastic _____ 21

Mattress, viscoelastic, short, with gynecological cutout _____ 21

Microsurgery wrist rest with integrated mounting bracket _____ 10

Microsurgery wrist rest (GENIUS EYE) _____ 11

Monitor holder, swiveling, on column _____ 26

Mounting adapter, 25 x 10 mm _____ 37

Mounting adapter, strengthened, 290 mm _____ 37

Mounting adapter with 90° redirection for leg holder _____ 17

Mounting bracket _____ 36

Mounting bracket SPECIAL “XL” _____ 36

Mounting bracket, swiveling _____ 12

Mounting bracket, swiveling (GENIUS) _____ 12

Mounting bracket, swiveling (VENUS) _____ 12

Mounting glider with clamping _____ 37

Mounting rail, 250 mm head section (GENIUS) _____ 38

Mounting rail, 340 mm foot section (GENIUS) _____ 38

Mounting rail, 250 mm, massiv seat section (GENIUS) _____ 38

Mounting rail, 340 mm, solid back and foot section (GENIUS) _____ 39

Mounting rail head section (GENIUS EYE) _____ 39

Mounting rail set (VENUS) _____ 39

P

Padded insert for arthroscopic surgery leg support _____ 19

Paper towel holder _____ 41

Paper towel holder, 550 mm (VENUS) _____ 41

Pedestal for foot control, type IV _____ 35

Protective film for splited foot section, visco (VARIUS) (1 pair) _____ 42

Protective film for foot section,visco (FOCUS, PRIMUS, GENIUS, VENUS) _____ 42

R

Radiolucent headrest, detachable (surcharge) _____ 9

Remote control holder, flexible _____ 33

Remote control holder on infusion stand _____ 32

Remote control holder, short _____ 33

S

Service kit _____ 43

Side extension/foot section extension, narrow (GENIUS) _____ 24

Side extension/foot section extension, narrow (OP) _____ 24

Side extension, radiolucent _____ 24

Side guard, hinged _____ 15

Side widener/foot section extension, narrow (GENIUS) _____ 24

Silicone head restraining belt for fastening on the headrest _____ 7

Standing plate (1 pair) _____ 19

Stomach positioning padding, trapezoid _____ 21

T

Tongue holder, extended, with flat plate rail _____ 11

Transport handles, foldaway, head end (1 pair) (GENIUS) _____ 40

Transport handles foldaway, head end, (1 pair) (FOCUS, PRIMUS, VARIUS) _____ 40

Transport handles foldaway, foot end (1 pair) _____ 40

Transport trolley for batteries _____ 45

Tray, Aluminium, for instrument table _____ 27

W

Wall mount for hand surgery table _____ 44

Wall mounting rail _____ 44

