

Dr's Light
AT

Best Choice for Best Treatment!

Dr's Light AT is a new conceptual curing light applied many advantages for user's demands.
Press the button and experience new light technology with Dr's Light AT.



-IC-WHT60
-IC-DCAU70
-IC-WHC50
-IC-WHW61
-IC-DCAU71
-IC-EZAU81

-CL-DLT30
-CL-DC20
-CL-DL10
-CL-AT24

-AL-DF10

'Best Choice for Best Treatment'

Dr's Light AT is a new conceptual curing light applied many advantages for user's demands.
Press the button and experience new light technology with Dr's Light AT.

Features

-   Focused light
4 curing light modes
-   Time setting available (Soft Start is fixed for 20 sec.)
330° head twist
-   Interchangeable battery pack
Big capacity battery



Without focus lens




With focus lens



Slim head



Disposable cover




Interchangeable battery

Specification


Model	CL-AT24
Dimensions (Included Battery)	Approx. 232 X 26 X 26 mm
Weight	Approx. 130 g
Wave length	440~490 (nm)
Power source	Lithium-ion Battery
Light source	High power LED

Accessory

Battery




Disposable Cover (approx. 100 pcs)



Comparison with Other products

Products (Dr's Light AT & Other products)



Various curing modes

Selectable light intensity and time as you want

